


# META FLEX

Percutaneous Flexible IM Nail System for Metacarpal and Metatarsal Fractures



# MetaFLEX

# Facilitating Improved Results for Both Metacarpal and Metatarsal Fractures

- **Simple, minimally invasive** procedure
- **Unique shaped design** to provide centering in varying diameters of IM canals
- **Color coded titanium wires** for easy identification of size
- **Cost effective treatment** of simple metacarpal or metatarsal fractures
- **Use of two MetaFLEX wires improves stability** over one slightly larger wire

## ORDERING INFORMATION


Cat. No.	Description
221.120	MetaFLEX Wire (1.2mm) x (100mm)
221.140	MetaFLEX Wire (1.4mm) x (100mm)
221.160	MetaFLEX Wire (1.6mm) x (100mm)
622.201	Complete Set




Metacarpal fixation




Metatarsal fixation




**Step 1:** Following a stab incision at the base of the metacarpal create an entry channel with a K-Wire/Awl in the metaphysis



**Step 2:** Insert the MetaFLEX IM Nail under fluoroscopic guidance advancing it into the medullary canal



**Step 3:** Reduce the fracture then continue to advance the MetaFLEX IM Nail under fluoroscopic guidance to the head of the metacarpal



**Step 4:** Cut the MetaFLEX IM Nail at the appropriate level



**SMALL BONE INNOVATIONS**

Small Bone Innovations, Inc.  
1380 South Pennsylvania Ave.  
Morrisville, PA 19067  
(215) 428-1791 / Fax (215) 428-1795  
SBI Customer Service: (800) 778-8837  
Technical Support: (866) SBI-TIPS  
www.totalsmallbone.com

Small Bone Innovations International  
ZA Les Bruyères - BP 28  
01960 Péronnas, France  
Tel: +33 (0) 474 21 58 19  
Fax: +33 (0) 474 21 43 12

Caution: United States federal law restricts this device to sale by or on the order of a physician.

7102 2M Copyright ©2006 Small Bone Innovations, Inc. All rights reserved.