


Frame Assembly and Application Just Got a Whole Lot Easier!



RingFIX™ Multi-hole Footplates and Tab Rings

- Lightweight components
- More options for easier application
- Strong constructs for Ilizarov Technique
 - Charcot Reconstruction
 - Arthrodesis
 - Ankle Distraction




RAISING EXPECTATIONS FOR LOWER LIMB SURGERY

RingFIX™

FOOT AND ANKLE FIXATION SYSTEM

More Options for Easier Application!

- Landmarks on the footplates and tabbed rings allow for quick and precise frame alignment.
- 3 rows of holes along the middle of the footplate allow for more options and easier placement of wires and half pins.
- Tab rings allow for easier modular application.



Landmarks allow for quick alignment – no hole counting required

Telescopic Rods available for easier distraction or compression.



More hole options where it counts!



Available in a wide variety of color coded sizes (140,160,180, and 210mm)



SMALL BONE INNOVATIONS

Caution: United States federal law restricts this device to sale by or on the order of a physician.

RingFIX is a trademark of Small Bone Innovations, Inc.
7007 xM Copyright ©2007 Small Bone Innovations, Inc.
All rights reserved.

Small Bone Innovations, Inc.
1380 South Pennsylvania Ave.
Morrisville, PA 19067
(215) 428-1791 / Fax (215) 428-1795
SBI Customer Service: (800) 778-8837
Technical Support: (866) SBI-TIPS
www.totalsmallbone.com

Small Bone Innovations International
ZA Les Bruyères - BP 28
01960 Péronnas, France
Tel: +33 (0) 474 21 58 19
Fax: +33 (0) 474 21 43 12

Anthony Viscogliosi
Founder and Chairman
AViscogliosi@totalsmallbone.com

MKT 20120 Rev C 9/07